






PURELL® Hand & Surface Sanitising Wipes



Antiseptic & disinfecting wipes for hands, surfaces & class I medical devices
Target markets: hospitals, dental clinics, health centres & emergency services

- Suitable for hands, surfaces and class I medical devices
- Proven to kill 99.99% of most harmful germs
- Skin friendly
- Designed and manufactured in France

TECHNICAL INFORMATION					
Composition	Ethanol CAS 64-17-5 (48 g/kg), Benzalkonium chloride CAS 85409-22-9 (3 g/kg), D-gluconic acid CAS 18472-51-0 (1,2 g/kg), fragrance				
Specifications	<ul style="list-style-type: none"> • Substrate: 100% Viscose (40 count), spunbonded (200 count, 225 count) • Fragrance: Menthol • pH: 4.5 – 6.5 				
Preservation	<ul style="list-style-type: none"> • Shelf-life: 30 months • Once opened use within one month • Store wipes in original packaging, with the container lid closed 				
Directions for use	1. Take a wipe	2. Wipe hands or surfaces for at least 1 min	3. Do not rinse, allow to dry	4. Close container after use	5. Single use wipe. Dispose of wipe in bin
					
Environment	<ul style="list-style-type: none"> • Dispose of the wipes in accordance with local regulations • Do not reuse the empty package 				
Regulation	<ul style="list-style-type: none"> • Biocidal product PT1 and PT2 • For professional use only 				
Precaution for use	<ul style="list-style-type: none"> • For the disinfection of hands, surfaces and medical devices, class 1 with smooth surfaces only • For more information refer to the product SDS • In case of emergency, contact the National Poisons Information Service • Use biocides safely. • Always read the label and product information before use. 				

ANTIMICROBIAL ACTIVITY OF IMPREGNATION SOLUTION

Active against	Standard	Contact time
Bacteria	EN 1040	1 min
	EN 1276* (v. March 2010)	30s
	EN 13727+A2*	30s
	EN 13697**	1 min
	EN 13697*	5 min
	EN 13727* : MRSA, <i>Enterococcus faecium</i> VRE, <i>Escherichia coli</i> ESBL	30s
	EN 14561*	5 min
	Yeast & Mould	EN 1275, <i>Candida albicans</i>
	EN 13624* (v.2013), <i>Candida albicans</i>	1 min
	EN 13697 *, <i>Candida albicans</i>	5 min
	EN 1650+A1, <i>Candida albicans</i>	1 min
	EN 1650, <i>Trichophyton mentagrophytes</i>	15 min
	EN 14562*, <i>Candida albicans</i>	5 min
Virus	Active against	
	H1N1 according to EN 14476+A1	30s
	Herpes simplex 1, Rotavirus	5 min
	PRV (surrogate for Hepatitis B virus) BVDV (surrogate for Hepatitis C virus), HIV	30s

*Study conducted in clean conditions **Study conducted in dirty condition

PACKAGING	PRODUCT CODE
40 x 40 count flowpack	92002-40-EEU
6 x 200 count canister	92106-06-EEU
6 x 225 count bucket	92206-06-EEU



GOJO Industries - Europe, Ltd.*
 Units 5 & 6
 Stratus Park
 Brinklow
 Milton Keynes
 MK10 0DE UK
 Tel: +44 (0) 1908 588444
 Fax: +44 (0) 1908 588445

Made in France by:
 Laboratoire Prodene Klint
 – a GOJO Family Company
 8 Rue Léon Jouhaux
 Z.A Pariest
 F-77435 MARNE-LA-VALLÉE – France
 Tél.: +33 (0)1 60 95 49 00
 Fax.: +33 (0)1 60 95 49 49