



INFORMATIVE NOTE

Name: REFLEX 99 NHR

Date: 06/04/2021 Rev: 016

MANUFACTURER

REFLEX SpA

Via Passeri, 2 46019 VIADANA (MN)

Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com

www.reflexx.com



Imballo - Pack

50		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R99 S	8032891632031	S	8032891637036
R99 M	8032891632048	M	8032891637043
R99 L	8032891632055	L	8032891637050
R99 XL	8032891632062	XL	8032891637067

DESCRIPTION: ■ disposable examination glove powder free, non sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff; long cuff ■ **Colour:** blue ■ **Feature:** internally chlorinated; textured surface and textured fingertips ■ **Packing:** 50pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES

Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE B - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 - EN 455 1,2,3 & 4

(CND T01020204 NON SURGICAL NITRILE GLOVES)



cat III

class I




Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/19. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length. min(mm)	Width (±5 mm)
R99/ S (7)	290	85
R99/ M (8)	290	95
R99/ L (9)	290	105
R99/ XL (10)	290	115
R99/ XXL (11)	290	125

Thickness mm	Palm (+/- 0.03)	0.15
--------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥6 N
------------------	------------------	------

EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass



EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015 TYPE B		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	9,6%
Hydrogen peroxide 30%	P	2	>30'	35,1%
Formaldehyde 37%	T	3	>60'	15,5%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD		Protection
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides
T Formaldehyde 37%	CAS: 50-00-0	Protection from the use of aldehydes

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In rare cases where transient hypersensitive skin reactions may occur, discontinue use and consult a doctor

■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment ■ Do not use gloves when there is a risk of getting caught in moving machine parts ■ **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com